STORZ

KARL STORZ — ENDOSKOPE

 20 5352 20-11x AUTOCON® II 400

 20 5352 20-11x AUTOCON® II 400

 20 5352 20-11x AUTOCON® II 400

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Instruction Manual

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Important information for users of KARL STORZ devices

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the AUTOCON® II 400. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

General warning

Users of KARL STORZ devices should be thoroughly trained in their use and the applicable medical procedures. Use of these devices and accessories should not be undertaken until the user has fully familiarized himself/herself with the instructions for use, assembly and care. Instruction manuals should be carefully studied and be available to the surgical team during the procedure; it is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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
Controls, displays, connectors and their uses

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Read the instructions carefully before operating
the equipment

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Connection socket for footswitch
(3-pedal footswitch 20013831 or
2-pedal footswitch 20013830 and
1-pedal footswitch 20013832)

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Service connection for software updates

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Potential equalization connector

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Applied part type CF as per IEC 601-1 (the neu-
tral electrode may remain applied to the patient
during defibrillation)

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* The SCB interface, based on the CAN field
bus, permits remote control of equipment func-
tions, as well as remote display of equipment
parameters.

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Dual/triple-pedal footswitch approved for cutting only (yellow footswitch)

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Dual/triple-pedal footswitch approved for coagulation only (blue footswitch)

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Yellow footswitch button activated for cutting and blue footswitch button activated for coagulation

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Single-pedal footswitch - only for coagulation

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+/- button for setting the cutting and coagulation parameters: performance and effect

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Indicator lights green when button is selected

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Return to level 2 (procedure overview)

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Other functions are assigned to this button

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For a detailed description of the symbols used for cutting and coagulation, see p. 40

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For a detailed description of NEs, see p. 22 and 58

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General warning	..II
Photographs of the equipment	..III
Controls, displays, connectors and their uses	..IV
Symbols employed	..V
General information	..2
The principle of high frequency surgery	..2
The AUTOCON® II 400 unit	..7
Safety instructions	..9
Warnings and cautions	..9
Intended use	..14
User qualification	..14
Safety precautions at the site of installation	..15
Safety precautions when operating the system	..15
Important notes on the use of the monopolar technique	..16
Safety devices	..21
Definitions and explanations for Autocon® II 400	..26
Definitions and explanations	..26
Symbols for effects which can be set	..40
Unit versions	..42
Installation and operating instructions	45
Unpacking the equipment	..45
Basic equipment	..45
Installing and connecting up the unit	..45
Operating the unit	..48
Maintenance	..69
Fuse replacement	..69
Cleaning and sterilization	..70
Maintenance	..75
Safety checks and inspections	..75
Servicing and repair	..76
Disposal	..76
Repair program	..77
Limitation of liability	..78
Manufacturer's warranty	..78
Technical description	..79
Troubleshooting	..79
Technical data	..82
Technical documentation	..104
General circuit diagram	..105
Spare parts, Recomm. accessories	..106
Spare parts list	..106
Recommended accessories	..107
Appendix	..108
Cleaning agents and disinfectants	..108
Electromagnetic Compatibility (EMC) Information	..110

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The principle of high frequency surgery


Biological tissue contains a more or less high concentration of electrolytes, making it sufficiently conductive to be treated electrosurgically. The thermal effect of HF current is used for separating (cutting) and coagulating tissue (desiccation of tissue and/or hemostasis).

HF currents must be used on the patient since low-frequency currents can stimulate nerve and muscle cells in the current flow due to electrochemical processes (electrolysis). These effects are small enough to be disregarded with frequencies above 100 kHz. In HF surgery, a distinction is made between the following application techniques:

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Monopolar operation techniques

With the monopolar operation technique, the electric current flows from the small-face active electrode to a large-face passive or neutral electrode, the patient's body always forming part of a closed circuit. Tissue separation or coagulation takes place at the active electrode.

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Bipolar operation techniques

With the bipolar operation technique, the current flows between two paired electrodes, the faces of which are either the same or different sizes, depending on the application (cutting or coagulation). Only a small part of the tissue is included in the circuit between the electrodes during application.

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Separation of tissue

When separating tissue (cutting), the cell sap of the tissue is vaporized. Due to the resulting vapor pressure, the tissue cells tear and the escaping vapor pressure causes the electrode to be insulated. If the voltage is high enough, the current flow is maintained by the commencement of arcing between the electrode and tissue. The thermal effect of the current causes the separation of tissue and the formation of surface slough, depending on the size of the arc.

Arcing commences at voltages of around 200 Vp. If voltages above 500 Vp are used, the result is substantial formation of slough due to the energy released in the arc.

For cutting, automatic control loops allow an optimal, reproducible cutting quality to be achieved, which is to the greatest possible extent independent of cutting depth, cutting speed, tissue resistance and such like. Power output is automatically adjusted to the level required for cutting.

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Tissue coagulation

To coagulate biological tissue, it must be heated to approx. 70°C.

To prevent the current from causing tissue separation, the power supply, and consequently the amount of heat, must be used specifically to produce denaturation of the tissue.

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The thermal effects which take place in biological tissue, depending on temperature, can be described as follows:

- Up to approx. 40°C: No visible cell damage
- Approx. 40°C and above: Reversible cell damage depending on the duration of current flow
- Approx. 49°C and above: Irreversible cell damage
- Approx. 70°C and above: Conversion of collagens into glucose – shrinking of tissue containing collagen – hemostasis
- Approx. 100°C and above: Dehydration/desiccation: phase transition of intracellular and extracellular water from liquid into gas. Adhesive effect of glucose after dehydration. Shrinking of the coagulum after dehydration.
- Approx. 200°C and above: Carbonization and fourth-degree pathological burns. Unpleasant odor of burnt tissue.

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- Standard Koagulation
- Forcierte Koagulation
- Spray-Koagulation

Due to the inhomogeneity of the electrical and thermal properties of tissue and the resulting uneven current density distribution within the tissue, the temperature rises with an irregular distribution when power is supplied.

The current density, and thus the temperature, is highest in the vicinity of the coagulation electrode. The temperature decreases very rapidly with increasing distance from the electrode so that coagulation is at first limited to the area surrounding the coagulation electrode. The effective contact area between the tissue and coagulation electrode has a decisive influence on the maximum achievable area of the coagulation zone.

Due to the spatial extent of the coagulation zone being dependent on the current density distribution, the technique with which the HF current is applied is also of importance. Here, a distinction is made between monopolar and bipolar surface application, monopolar and bipolar puncture application as well as monopolar and bipolar grip application.

Coagulation modes

The spatial extent of the coagulation zones depends on the type of energy supply and application technique.

It is useful to differentiate between 3 coagulation modes:

- Standard coagulation
- Forced coagulation
- Spray coagulation

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

Standard coagulation

With standard coagulation, the generator voltage is kept so low that no arcing is able to take place between the electrode and tissue during the coagulation process. Standard coagulation is recommended for all instances of coagulation in which monopolar or bipolar coagulation electrodes are used which are brought into direct contact with the tissue that is to be coagulated. Standard coagulation provides for deep coagulation while avoiding carbonization of the tissue.

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Forced coagulation


With forced coagulation, electric arcs are intentionally permitted between the coagulation electrode and tissue to enable the formation of surface sloughing with nevertheless deep reaching coagulation. This particularly applies if only small-face electrodes are available for coagulation. To avoid the risk of a cutting effect, a modulated HF voltage must be used. Tissue sloughing will, however, occur.

A typical application is TUR, in which the cutting electrodes are also used for coagulation purposes.

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
Spray coagulation

If sufficiently high voltages pulses are used for coagulation, arcs occur between the electrode and tissue of such a length that contact of the tissue with the electrode is no longer necessary. Spray coagulation is primarily used for surface coagulation or for hemostasis of vessels which cannot be directly reached with a coagulation electrode (e.g. in bone fissures).

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Legend:

- A: Power required for cutting
 B: Power output without control of unit
 C: Excess power

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The operator bases his assessment of the quality of the tissue cutting on the various parameters of the entire system. Such parameters include the geometry and size of the instrument used, the nature of the tissue at the site of the operation, the penetration depth and cutting speed of the instrument, the conductivity of the irrigant used, and the application of the neutral electrode etc.

These 'external' parameters, which can be extremely variable, influence the properties of conventional HF equipment to such a degree that the operator must constantly readjust the power or generally set a higher power level. The excess power, however, produces an unnecessary degree of surface slough formation and exposes the patient to unnecessarily high current and the associated risks.

Fig.: Power output (P) in relation to time (t) when cutting tissue.

The cutting action during tissue cutting is influenced to a great extent by the electric arc maintained between the instrument and the tissue. This realization led to the development of a new generation of high frequency surgical units.

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The AUTOCON II 400 unit

The AUTOCON® II 400 HF surgical unit is a universally applicable unit which permits the precise setting of various cutting and coagulation qualities as well as reliable reproducibility of these settings. Due to automatic open-loop and closed-loop control of the HF voltage, the HF voltage and intensity required for achieving the desired cutting quality are provided in each phase of the incision.

The AUTOCON® II 400 permits both monopolar and bipolar incisions and provides a suitable cutting and coagulation mode for every application.

The bipolar coagulation mode additionally permits operation with the autostart and autostop functions.

In conjunction with automatic open-loop and closed-loop control of cutting and coagulation quality, the various safety devices of the instrument provide an exceptionally high level of user and patient safety. Software-supported test programs ensure easy and rapid servicing.

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Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be read carefully to ensure the safe and effective operation of this product. To make the words Warning and Caution stand out more clearly, they are accompanied by a pictogram.

Definitions

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

Note: A Note indicates special information about operating the unit.

Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

Warning: The electrical installation of the operating room must comply with the applicable IEC requirements.

Warning: Keep out of reach of patients.

Warning: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

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Warning: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

Warning: Grounding reliability can only be achieved when the equipment is connected to a "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly. The unit should not be used if any damage is evident.

Warning: Risk of explosion if used in the presence of flammable anesthetics.

Warning: Always unplug the unit before carrying out any maintenance work on it (e.g. cleaning).

Warning: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

Warning: Contact with noninsulated instruments by the active high frequency electrode can lead to unintentional burns. We therefore recommend using insulated instruments in the vicinity of the active electrode.

Warning: Avoid currents flowing via tissue bridges. The resulting current concentrations at the base adjoining healthy tissue may cause burns.

Warning: Use an insulated eyepiece on the endoscope to reduce the risk of capacitive leakage currents flowing to the user or to the attached camera system.

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Warning: Please always observe the maximum permitted high frequency peak voltage for operation of the instruments. Failure to observe the maximum permissible high frequency peak voltage for the operation of accessories to high frequency generators can lead to destruction of the insulation and thus to possible burns at the contact points of the instrument to the patient or user.

Warning: Disturbances produced during operation of the high frequency surgical device can impair the function of other electrical devices.

Warning: The unit may only be operated under the environmental conditions specified.

Warning: Please ensure that there is sufficient air circulation for your unit.

Warning: Please always observe cleaning and sterilization instructions.

Warning: The unit must not be operated if the touch screen or display is defective

Warning: It is not permitted to press more than one touch point simultaneously.

Warning: During simultaneous use of high frequency surgery and physiological monitoring equipment on a patient, every monitoring electrode should be positioned as far away as possible from the active and neutral electrodes of the high frequency surgical device. **Die Verwendung von Nadelelektroden für die Überwachung ist zu vermeiden.**

In all cases monitoring electrodes which contain devices for limiting the high frequency current are recommended.

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Use of several HF units on one patient

Warning: Only use units with floating F application parts. Always operate a unit with its own neutral electrode.

Avoid high frequency electric circuits, i.e. ensure adequate spacing (i.e. at the top, bottom, in the centre, no switching of active to neutral electrodes).

All units should be operated according to the instructions provided in the respective instruction manuals.

Injury due to incorrect high frequency parameter

Warning: Always select a high frequency output power setting suitable for the particular intended procedure in order to avoid, on the one hand, thermal invasion of the tissue, which can be caused by too high a setting, or, on the other hand, insufficient coagulation resulting in excessive bleeding, which can be caused by too low a setting.

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Safety instructions**Warnings and cautions**

Warning: Avoid explosive gas concentrations that may be present in the area where a high frequency surgical procedure is being conducted. The following situations, among others, must be taken into consideration:

Gas which may support combustion is sometimes present in the gastrointestinal tract of an unprepared patient. Certain substances used for the preparation of the lower GI endoscopy can enhance methane production. This is particularly relevant in colonoscopy, but has also been recorded in the upper GI tract. In addition, during transurethral resection of the prostate, it has been recorded that hydrogen can accumulate in the bladder above the irrigation solution, this is particularly true when a strong arc is used.

Working with explosive gas mixtures may lead to deflagration/explosion.

Caution: Only operate the unit with the voltage stated on the unit's identification plate.

Caution: Use only fuses of the correct rating.

Caution: Avoid allowing liquids to enter the unit. Do not store liquids on or above the unit. It must also be ensured that no liquid bottles or tubing are placed near the unit.

Caution: The signal input and signal output parts of this equipment are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

Caution: It is essential to refer to the Instruction Manual for proper disinfection/sterilization instructions.

Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Note: In the case of HF generators, attention must be paid to the maximum peak voltage allowed. Please observe the instruction leaflet provided to prevent damage to the material.

Note: Inspect the electrode leads and endoscopic accessories regularly for damage to the insulation.

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Before using the unit on the patient for the first time it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The AUTOCON® II 400 is a high frequency surgical unit for monopolar and bipolar cutting and coagulation.

The AUTOCON® II 400 may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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User qualification

The AUTOCON® II 400 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the instrument. Such instruction should only be provided by persons equipped to do so on the basis of their skills.

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Patient profile

The use of Autocon II 400 is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

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Affected body part or tissue type

The field of application is not confined to certain parts of the body or types of tissue.

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User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Adequate knowledge of the contents of the instruction manual
- Adequate knowledge of the purpose and whereabouts of the medical products logbook
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

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Use

The device is used exclusively in hospitals and physicians' practices under the ambient conditions given in the technical data (page 98).

Other conditions of use:


- Frequency of use up to several times a day
- Length of use from a few minutes to several hours a day
- Position it on a vibration-free, fixed base
- Can be moved only by installing it on an equipment trolley
- Can be swiveled only by installing it on ceiling supply units
- **Intervallbetrieb wie auf dem Typenschild angegeben. 25% relative Einschaltdauer (z.B. 10 sec aktiviert, 30 sec. deaktiviert)**

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Safety precautions at the site of installation

The unit may only be used in medical rooms the installations of which conform to applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures of these, the unit must not be operated inside the demarcated hazard zone. Examples of such substances are: anesthetic ether (diethyl ether, cyclopropane) as well as combustible, volatile skin cleansers and skin disinfectants which may create an explosive atmosphere (e.g. detergent ether, petroleum ether). The unit is equipped with a connector for attaching a ground line. It should be connected in accordance with national regulations.

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Safety precautions when operating the system

It is the user's responsibility to make sure the unit is safe and operates properly before using it.

During treatment with the AUTOCON® II 400 the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if visual observation of the action of the unit is ensured.

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Important notes on the use of the monopolar technique

Risk of burns

HF surgery always entails certain risks for both user and patient. Above all, contact burns can occur due to unintentional thermal damage of tissue. With a knowledge of the possible causes, and by observing the safety instructions below, burns can, however, be avoided.

Burns due to incorrect application of the neutral electrode

With incorrect or no application of the neutral electrode, there is a high risk of thermal tissue damage both at the point of application of the neutral electrode and at other points of the patient's body. When applying neutral electrodes, it is therefore imperative that care be taken to ensure that the effective contact area or conductance between the neutral electrode and the patient is adequate for the HF power used or the intensity of the HF current. Here, effective contact area refers to the area of the neutral electrode in electroconductive contact with the patient's skin during HF surgery. When selecting the neutral electrode, take the following empirical values into account:

Neutral electrode made of electroconductive silicone:	1.5 watts/cm ² or 10 mA _{rms} /cm ²
Neutral electrode made of metal without conductive jelly:	3.0 watts/cm ² or 15 mA _{rms} /cm ²
Disposable neutral electrode:	5.0 watts/cm ² or 20 mA _{rms} /cm ²

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The neutral electrode is secured to either the upper arm or thigh, depending on which is nearest the operating field. To ensure that the entire face of the neutral electrode is applied, it must be securely fastened with rubber straps and then wrapped with a bandage so that it cannot slip, even if the patient should move (note the risk of decubital necrosis). The neutral electrode must always be in perfect condition. Grease, oxidation, dirt etc. on the contact surface reduces the size of the electroconductive area and can lead to dangerously high current densities.

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Warning: A short circuit fault is detected when using a double-faced NE. Power is disconnected until optimal conditions return. Subsequently, a renewed HF activation is essential.

Warning: The attachment of a neutral electrode made of conductive silicone (27805) prevents monitoring.

Note: The HF current flows to the neutral electrode through the tissue. It should be remembered that the distribution of current at the proximal corners and edges to which the current flows may be higher than at the distal corners and edges. For this reason, when arranging the neutral electrode care must be taken to ensure that the longer side of the neutral electrode is facing the operating site.

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
Patients must be correctly positioned before performing any operations using HF surgical instruments. An electrically insulating, thick and dry pad must be placed between the patient, operating table and clamps. If moisture, perspiration etc. are anticipated during the operation, penetration of moisture through this pad must be prevented by waterproof sheeting or other suitable measures. The patient's body tissue must not be allowed to come into contact with nearby electrically conductive objects (examples: operating table, metal clamps, wet towels, etc.).

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Since contact areas will generally be small in such cases, burn-necroses could arise from high frequency leakage currents due to high local current densities. In cases where contact between the neutral electrode and patient is poor this effect is increased.

Areas subject to more severe perspiration, extremities resting against the body or skin-to-skin contact should be dried by placing cover cloths in-between. Urine should be drained off by a catheter. If an injection cannula is used as an ECG electrode, the metal cone must not rest on the skin; the same applies to the leads of monitoring equipment.

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Burns due to unintentional activation of the HF generator

Patients can also suffer burns due to unintentional activation of the HF generator. Unintentional activation can occur, for example, if the footswitch is inadvertently depressed, if there is a short circuit in one of the cables, if electrically conductive liquids enter the electrode handle or the footswitch, or if the HF generator is defective. To prevent such occurrences:

Warning: To prevent such occurrences:

- never leave electrodes lying on or near the patient in such a way that could result in direct or indirect (through the mediation of electrically conductive objects) contact between electrode and patient,
- ensure that the acoustic warning signal indicating the operational status of the HF generator is always set loud enough to be easily heard.
- ensure that the visual warning signal can be seen by all attending personnel.

Warning: Electrodes which are not currently in use must not come into contact with the patient.

Warning: Skin-to-skin contact should be avoided.

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Burns from contact with hot electrodes

Cutting and coagulation electrodes can become very hot during cutting or coagulation procedures. Burns can result from active or inactive electrodes coming into contact with tissue immediately following cutting or coagulation procedures.

Burns due to improper setting of HF power output

The optimal HF power level setting is that which yields the best results in terms of the desired effect (cutting/coagulation) with minimal incidence of undesirable side effects (i.e. carbonization of tissue or sparking). Since undesirable side effects generally occur at excessively high power levels and/or if effects are too high, HF power levels should always be kept as low as possible.

If the desired effects should suddenly no longer be attainable at power levels which experience has shown to be optimal, always check to be certain that everything is in proper working order before increasing the power level and effects:

- Is the neutral electrode properly connected?
- Is the active electrode properly connected?
- Are active and neutral electrodes clean?
- Is the HF generator power level and/or effect correctly set?

Burns arising from ignition of combustible media

Skin cleansing agents, disinfectants, and anesthetic gases, as well as endogenic gases, may be readily combustible and can even explode. Since sparking inevitably occurs in HF surgery, vapors or gases evolved by such media could be ignited and result in burn injuries to patients. Thus, before switching on the HF generator, make certain that all combustible media and all objects that have been in contact with them (tissues, swabs, etc.) have been removed from the operating area.

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Burns due to unsuitable application

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Stimulation of nerves and muscles

Due to low-frequency electric currents, nerves and muscles of the patient can be stimulated during HF surgery. The arcs between the electrode and tissue, which are unavoidable for cutting processes, spray coagulation and standard coagulation, have the effect of rectifying part of the HF alternating current, whereby low-frequency current components occur which are able to cause nerve stimulation. When using HF surgery on electrically stimuable structures, contractions of the muscles concerned must therefore be expected.

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Heart pacemakers

Irreparable damage to heart pacemakers and interference with pacemaker operation, which could lead to ventricular fibrillation, are likely to result if high frequency surgery methods are used on patients with implanted heart pacemakers or pacemaker electrodes.

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Implants

In the case of patients with electrically conductive metal implants it must be ensured that with all applications involving HF surgery the current flow or current path of the HF electric current does not pass through or over the implant. At the current entry and exit points into and from the bio-tissue to the implant, current concentration may cause irreparable internal burns.

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Safety devices**Safety circuit in case of output error caused by the unit**

Deviations from the setpoint of the respective HF output that has been set are corrected by the unit. An independent safety circuit monitors the power output and switches it off if proper cutting or coagulation is no longer ensured.

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Monitoring the activation time

With patient application, the HF generator is only activated briefly when used as prescribed. The activation time of the HF generator is therefore monitored. If the maximum activation time is exceeded, this is indicated by the error message A-1C on the display. The power output is then switched off. This safety circuit reduces damage resulting from unintended activation of the electrodes. Limitation of the maximum activation time is 99 seconds on delivery of the unit. Specific selection can be made under 'Application Settings' ranging from 0 to 99 seconds and unlimited.

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Note: Continuous activation can cause irreparable damage to the unit.

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Safety Devices

Monitoring the neutral electrode

The AUTOCON® II 400 has a safety circuit which monitors the electrical connection between the neutral electrode and unit and, when using double-faced electrodes, the correct application on the patient. The neutral electrode (NE) status display indicates how well the NE has been positioned on the patient. If it is positioned correctly on the patient a green neutral electrode is displayed. If a critical limit value is too low, this is displayed in red.

Warning: If no neutral electrode is connected, the symbol for neutral electrode lights up in red. If this is not the case, the unit must not be used!

If a **single-faced** neutral electrode is used, the status display lights up green when the unit and the neutral electrode are correctly connected. **Bei einer einflächigen NE wird die Kontaktqualität überwacht (NE Widerstand), nicht jedoch die Stromdichte oder Unsymmetrie.** Partial detachment of a single-faced neutral electrode always increases the current concentration under the neutral electrode and results in thermal damage to tissue due to the limited contact area.

When using **double-faced** neutral electrodes, the electrical resistance the current density and the symmetry between the two contact areas of the neutral electrode is measured. If the resistance is outside of a tolerance range, this is recognised by the unit and a signal is given.

If the connection is interrupted or the neutral electrode positioned poorly (NE symbol lights up in red) the power output in the monopolar operating modes is blocked.

If an attempt is made to activate the HF generator in a monopolar operating mode, an acoustic warning signal is given. At the same time, an error message is displayed on the screen.

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Self-test program

Each time the unit is powered up, a self-test program is run. First the logo is displayed for approx. 3 seconds, afterwards a white screen appears for approx. 2 seconds. This serves to optically assess if the screen is functioning correctly. If any part of the white page shown displays another color, the screen is defective and the unit should no longer be used. This is for safety reasons as the HF output data may might not be displayed or displayed incorrectly. Please notify your KARL STORZ service partner at once!

If the unit detects a fault, a visual and, **abhängig von der Klassifizierung, auch acoustic** warning is given. The error codes can be found in the table 'AUTOCON® II 400 Error messages' on page 80. After the self-test, a procedure overview appears on the monitor. The date and time of day are displayed on the bottom right. The flashing seconds display on-screen shows that the function is activated (alive signal)

Switch-off option for HF activation

HF activation can be switched off for every HF socket and every mode by selecting 'Cut off' or 'Coag off' Activation using the footpedal and/or fingerswitch is then no longer possible.

This is a safety function which prevents a cutting mode from being accidentally activated for example during surgery where only coagulation is required.

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Safety instructions Safety Devices

Alarm system

The AUTOCON® II 400 is equipped with an alarm system. The alarm messages are displayed on the display in color and backed up with an acoustic signal.

The visual and acoustic messages used are broken down into three priority categories:

Alarm category	Color	Flashing frequency
High priority	Red	2 Hz/0,5 sec
Medium priority	Yellow-orange	0.5 Hz/2 sec
Low priority	Turquoise	Constantly on

Depending on the priority the alarm message is backed up with an acoustic signal:

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Low priority: No acoustic backup signal.

Note: Any low priority signal will be overwritten by an alarm signal with a higher priority.

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
The alarm system is capable of identifying the alarm with the highest priority in the status bar. In the event of more than one alarm with the same priority, only the alarm which occurred first is shown in the status bar. The alarm system classifies the alarms into those which do not permit further use of the unit and those which permit further usage.

In addition to alerts, there are also purely informative messages. These messages always appear on the user interface for just 3 seconds.

The list of all error messages relevant for the user and corresponding priorities can be found in the table "AUTOCON 400 Error Messages" on page 80. Errors which can only be rectified by technical service personnel are listed in the service manual.

Note: The user's position must always be selected so that any alarm message can be immediately identified. The distance from the unit should not be more than 4m.

Note: The sound pressure of the audible alarms is at least 65 dBA at a 1 m distance from the front of the unit in the horizontal direction.

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Sockets

The unit can be equipped with four different HF output sockets. At the factory, several unit versions are provided with monopolar, bipolar, multi-function and neutral electrode sockets. Details on this can be found here in the table under item 'Unit versions' (p. 43). With the exception of the neutral electrode socket, information about the general monopolar, bipolar or bipolar multi-function types and the equipment location is available. The user can assign a socket with the following parameters:

- cutting mode
- coagulation mode
- cutting and coagulation activation type

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Modes

A mode is characterized by different generator settings such as frequency, peak voltage, modulation, etc. A total of 10 different cutting and 9 different coagulation modes are available. **Nicht jeder Modus ist auf jeder Buchse abrufbar.** The user can also assign the following parameters to a mode:

- power limitation cutting,
- cutting effect,
- power limitation coagulation
- coagulation effect.

If a new mode is selected, this will start with specific setpoints for power limitation and effect. The modes 'Cut off' and 'Coag off' represent a special case here. Activation of these modes blocks activation via the footswitch or fingerswitch.

The default settings of the operating mode are described on page 50.

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Power limitation

The HF power output limit is max. 370 W +8% / -20%. Maximum limit that can be set by the user is 300 W. However, this maximum limit is not available in every mode.

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Effect

The cutting or coagulation effect can be preselected by the user in several intensities. The number of intensity levels depends on the mode selected. (see table p. 41)

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Activation types

Three footswitch types are available to activate the HF energy. Single-pedal, dual-pedal and triple-pedal footswitches. Activation is assigned to an output socket. For monopolar sockets there is the option of using the fingerswitch of certain instruments for activation. This function is always activated and can be switched off with the modes "CUT OFF" and "COAG OFF". In certain bipolar modes the activation types "AUTOSTART 1" and "AUTOSTART 2" are also available. When contact is made with the tissue, depending on the impedance, the HF generator is automatically activated after a time delay defined by the user.

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Buchse	Schneid-Modi	Coag-Modi
Bipolar	Bipolar Cut	Bipolar Soft Coag
	Cut off	Bipolar Soft Auto Stop Coag aus
Bipolar Multifunktion	Saline-C-Cut	Bivascular Safe
	Saline-Time-C-Cut	Saline Coag
	Cut aus	Saline-Time-Coag Coag aus
Monopolar	Power Cut	Standard Coag
	Top Cut	Forced Coag
	Lap-C-Cut	Spray Coag
	C-Cut	Coag aus
	Gastro-Cut	
	Papillo Cut	
	Cut aus	
NE	Widerstandsmessung	

Socket	Cutting modes	Coag modes
Bipolar	Bipolar Cut	Bipolar Soft Coag
	Cut off	Bipolar Soft Auto Stop Coag off
Bipolar Multifunction	Saline-C-Cut	Bivascular Safe
	Saline-Time-C-Cut	Saline Coag
	Cut off	Saline-Time-Coag Coag off
Monopolar	Power Cut	Standard Coag
	Top Cut	Forced Coag
	LAP-C-Cut	Spray Coag
	C-Cut	Coag off
	Gastro-Cut	
	Papillo Cut	
	Cut off	
NE	Resistance measurement	

Разъем	Режим резания	Режим коаг.
Биполярный	Bipolar Cut	Bipolar Soft Coag
	Cut off	Bipolar Soft Auto Stop
Биполярный многофункциональный	Saline-C-Cut	Bivascular Safe
	Saline-Time-C-Cut	Saline Coag
	Cut off	Saline-Time-Coag Coag off
Монопольный	Power Cut	Standard Coag
	Top Cut	Forced Coag
	LAP-C-Cut	Spray Coag
	C-Cut	Coag off
	Gastro-Cut	
	Papillo Cut	
	Cut off	
NE	Измерение сопротивления	

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Warning: The mode and power setting should be adapted to the intervention and should be as low as possible.

Warning: When performing monopolar HF surgery on bio-tissue with fine structures, there is the risk of burns due to local current concentration. For this reason, bipolar HF surgery is recommended in these areas.

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POWER Cut

Properties

Reproducible voltage-controlled atraumatic cutting, minimizing trauma to tissue, low to medium level of hemostasis.

The POWER Cut mode is equipped with PPS (Power Peak System). The AUTOCON® II 400 is equipped with an automatic power control system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arc necessary for the cutting quality selected even in the case of low-resistance loads. This is particularly important when the cutting wire electrode is pressed firmly against the tissue to be cut before activating the HF generator so that the cutting wire electrode covers a relatively large area and is therefore in low-resistance contact with the tissue. This is, for example, the rule with TUR and endoscopic polypectomy. In such cases the HF generator has to provide an above-average level of power on activation to ensure incision can be made without delay. Thanks to this device the average power can also be limited to relatively low levels, thus reducing the risk of unintentional thermal damage to tissue.

Fields of application:

All cutting procedures in tissue with good electric conductivity: e.g. muscle tissue, vascularized tissue. Dissection or cutting fine structures.

Suitable electrodes:

Needle electrodes, knife electrodes, spatula electrodes, cutting loops.

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TOP Cut

Properties

Arc-controlled reproducible cutting which is atraumatic to tissue, in particular for tissue varying in nature and with poor electric conductivity.

PPS (Power Peak System)

The TOP Cut mode is equipped with PPS. The AUTOCON® II 400 is equipped with an automatic power control system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arc necessary for the cutting quality selected even in the case of low-resistance loads. This is particularly important when the cutting wire electrode is pressed firmly against the tissue to be cut before activating the HF generator so that the cutting wire electrode covers a relatively large area and is therefore in low-resistance contact with the tissue. This is, for example, the rule with TUR and endoscopic polypectomy. In such cases the HF generator has to provide an above-average level of power on activation to ensure incision can be made without delay. Thanks to this device the average power can be limited to relatively low levels, thus reducing the risk of unintentional thermal damage to tissue.

Fields of application

Includes the cutting of fatty structures and cutting under water, e.g. with TUR-P.

Suitable electrodes

Knife electrodes, spatula electrodes, cutting loops.

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LAP-C-Cut

Properties

High level of hemostasis with slightly arrested cutting.

Fields of application

E.g. incisions during laparoscopy or endoscopic interventions which require excellent primary hemostasis during cutting and tolerate slightly arrested cutting.

Suitable electrodes

Electrodes with a large contact surface, knife and spatula electrodes, hook electrodes.

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C-Cut

Properties

High level of hemostasis with slightly arrested cutting.

Fields of application

E.g. cutting with low levels of bleeding during monopolar endoscopy with TUR-P and TUR-B, hysteroscopy.

Suitable electrodes

Electrodes with a large contact surface: knife and spatula electrodes, cutting loops and ribbon snare electrodes.

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PAPILLO Cut

Properties

Cutting consists of alternating cutting and coagulation phases. Cutting is easily controlled and is characterized by a reproducible and preselectable coagulation property that is specially suited to papillotomy during cutting.

Fields of application

Endoscopic interventions in which alternating cutting and coagulation is required with a single activation.

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Cutting current which is used with conductive NaCl irrigation solution in combination with bipolar KARL STORZ resectoscopes for cutting with low levels of bleeding during TUR-P and TUR-B and also surgical hysteroscopy with conductive NaCl irrigation solution

PPS (Power Peak System)

Fields of application

E.g. cutting with low levels of bleeding during TUR-P and TUR-B and surgical hysteroscopy.

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During use, the Saline-Time-C-Cut is emitted with automatic timed switch-off from 0.1 to 0.9 seconds. With premature disabling of the footpedal, the power output is immediately interrupted. During the activation time, the Saline-Time-C-Cut current is emitted when the yellow footpedal is depressed.

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Definitions and explanations for AUTOCON® II 400

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BIPOLAR Cut

Properties

Cutting current which only flows directly around the distal end of the applicator. The effect intensities can be used to vary the degree of hemostasis at the incision edge.

PPS (Power Peak System)

The BIPOLAR CUT mode is equipped with PPS. The AUTOCON® II 400 is equipped with an automatic power control system that detects low-resistance loads and controls the HF generator so that it briefly provides sufficient power to ensure the HF voltage or intensity of the electric arc necessary for the cutting quality selected even in the case of low-resistance loads. This is particularly important when the cutting wire electrode is pressed firmly against the tissue to be cut before activating the HF generator so that the cutting wire electrode covers a relatively large area and is therefore in low-resistance contact with the tissue. This is, for example, the rule with TUR and endoscopic polypectomy. In such cases the HF generator has to provide an above-average level of power on activation to ensure incision can be made without delay. Thanks to this device the average power can be limited to relatively low levels, thus reducing the risk of unintentional thermal damage to tissue.

Suitable electrodes

Special applicators (bipolar electrodes with a rigid or height-adjustable cutting needle) during laparoscopy, neurosurgery and ENT.

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STANDARD Coag

Properties

High level of slough formation in tissue is prevented, adhesion of electrode to tissue is reduced. If you wish to make use of the potentially great 'coagulation depth' of SOFT Coag, select a low effect intensity and coagulate for an extended period of time. If you can only coagulate for a brief period, select a high effect intensity.

Fields of application

For almost all operations which require fast and reliable coagulation or where adhesion of the electrode would have an adverse effect on the coagulation process.

Suitable electrodes

Electrodes with a large contact area, e.g. ball electrodes for deep coagulation.

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Fields of application

Contact coagulation, or clamp coagulation via insulated monopolar forceps for example.

Suitable electrodes

Ball electrodes for contact coagulation. Insulated monopolar forceps for clamp coagulation. Lancet electrodes via lateral contact with a large contact surface.

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BIPOLAR SOFT Coag

Properties

Low voltages, slough formation of tissue is prevented, adhesion of electrode to tissue is greatly reduced. If you wish to make full use of the potentially great coagulation depth of BIPOLAR SOFT Coag, select a low effect intensity and coagulate for an extended period of time. If you can only coagulate for a brief period, select a high effect intensity. In comparison with other Coag modes you will still achieve a great coagulation depth, although you will not make full use of the potential coagulation depth of BIPOLAR SOFT Coag.

BIPOLAR SOFT Coag with AUTOSTART

In the window Activation type you can select the AUTOSTART function for BIPOLAR SOFT Coag. This additional function starts coagulation automatically after a defined time period when the instrument touches the tissue. Touching conductive metal material with the bare tips of the forceps, for example, will always initiate activation of the HF current with the additional AUTO START function!

Warning: When introducing bipolar instruments into a metal trocar, the Auto Start function should be disabled beforehand!

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BIPOLAR SOFT Coag with AUTOSTOP

Properties

Identical with BIPOLAR SOFT Coag; however, it is firmly coupled with the AUTOSTOP function.

The additional AUTOSTOP function ends activation automatically before the tissue adheres to the instrument or slough formation commences.

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SALINE-Coag

Properties

Coagulation mode for use with conductive NaCl irrigation solution in combination with KARL STORZ bipolar resectoscopes for hemostasis with TUR-P and TUR-B interventions as well as hysteroscopy. Here, the active resection loop must make contact with the bleeding vessel and the HF current activated at the same time.

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SALINE-TIME-Coag

Properties

Used for timed output of Saline-Coag HF current from 0.1 - 0.9 sec. Can also be interrupted prematurely by the user by being disabled.

Fields of application

E. g. for hemostasis with TUR-P and TUR-B interventions as well as hysteroscopy.

Suitable electrodes

Cutting loops recommended by KARL STORZ, vaporization electrodes, needle electrodes for bipolar resectoscopes from KARL STORZ with NaCl irrigation solution.

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Bi-VASCULAR-SAFE

Properties

The mode which is available for the unit versions -113 to -116 at the MF socket (multifunction socket) is designed for the coagulation of vessels up to approx. 4mm in diameter. AutoStop deactivates automatically before the tissue can adhere to the instrument. A brief audio signal is heard when it is disabled.

Fields of application

Open surgery, laparoscopy and endoscopy.

Suitable electrodes

Special bipolar KARL STORZ instruments with an MF connecting cable.

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Operating mode	Effect intensity								Image in socket information display	Image in selection button
	1	2	3	4	5	6	7	8		
Cut/Coag off										
TOP-Cut										
POWER-Cut										
Bipolar-Cut										
C-Cut										
Lap-C-Cut										
Saline-C-Cut + Saline-Time-C-Cut										
Gastro-Cut										
Papillo-Cut										
Standard Coag										
Forced Coag										
Spray Coag										
Bipolar Soft Coag										
BI-Soft Autostop										
BI-Vascular-Save										
Saline-Coag + Saline-Time-Coag										

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





Unit versions:

For the AUTOCON® II 400 the following equipment versions are available with different socket

AUTOCON® II 400 Set	AUTOCON® II 400	Socket Position			
		1	2	3	4
Standard					
20 5352 01-111	20 5352 20-111	Bipolar 	Bipolar 	Unipolar 3-pin and Erbe 	NE 6.3 mm jack 
20 5352 01-112	20 5352 20-112	Bipolar US-2-pin 	Bipolar US-2-pin 	Unipolar 3-pin and Bovie 	NE 2-pin 
Bipolar					
20 5352 01-113	20 5352 20-113	Bipolar 	Bipolar 	Bipolar Multifunction 	Blind Socket 
20 5352 01-114	20 5352 20-114	Bipolar US-2-pin 	Bipolar US-2-pin 	Bipolar Multifunction 	Blind Socket 
High-End					
20 5352 01-115	20 5352 20-115	Bipolar 	Bipolar Multifunction 	Unipolar 3-pin and Erbe 	NE 6.3 mm jack 
20 5352 01-116	20 5352 20-116	Bipolar US-2-pin 	Bipolar Multifunction 	Unipolar 3-pin and Bovie 	NE 2-pin 

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Unpacking the equipment

Carefully remove the AUTOCON® II 400 and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Installing and connecting up the unit

Warning: The AUTOCON® II 400 HF surgery unit including the accessories connected may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes. This unit is not intended for use in hazardous zones. Do not operate the unit within the demarcated hazard zone while explosive anesthetic gases are in use.

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Place the unit on a level surface.

Warning: HF surgery equipment generates radio-frequency and electromagnetic fields that can interfere with the operation of sensitive electronic equipment. HF surgical equipment, and particularly its cabling, should thus never be installed close to any equipment sensitive to electromagnetic interference.

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
The instrument is equipped with a connector ⑧ for attaching a potential equilization. The instrument's ground line should be installed by a qualified electrician.

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Checking selected line voltage (115/230 V)

Note: Remove the power cable before opening the fuse holder.

The Autocon® II 400 units are available in two voltages: 115 V and 230 V. You can tell the 115 V variants by the 'U' in the article number and they may only be operated at 115 V. The 230 V variants have a '-' in the article number on the manufacturer's identification plate and may only be operated at 230 V. How to check whether the fuse ratings and the voltage selection module are correct is described below.

If necessary, remove the cover of the line fuse holder/voltage selection module (see 'Fuse replacement' section). Take out the line fuse holder/voltage selection module. Turn through 180° until the indicator (115 or 230 V) matches the local voltage supply.

The setting 115 V is valid for a voltage range of 100-120 24 VAC, and the setting 230 V for a range of 230-240 VAC.

Make sure that the fuses used are suitable for the line voltage selected and are as stated on the manufacturer's identification plate:

20535220-11x, 230 V: T4 AH250V
(Art. no. 2027690)

20535220U11x, 115 V: T8 AH250V
(Art. no. 2028090)

Insert line fuse holder/voltage selection module back into unit, close cover.

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Connect power cord when the unit is switched off. Insert power cord into power cord receptacle ⑨ as far as it will go.

Warning: Make sure to use the correct line voltage and power cord type.

Note: The instrument may only be operated with the power cord delivered by the manufacturer or a power cord of similar quality ('Hospital Grade Plug'), which has a national inspection seal.

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Only operate the unit with the voltage stated on the identification plate.


Warning: Only connect the power plug to the power supply away from areas where there is danger of explosion.

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
Insert the connecting cable of the footswitch into socket ⑤. Secure cable using screwed connection (union nut).

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SCB

Note: To prevent the SCB connecting cable being pulled out accidentally, the SCB connector possesses a protection device. Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑦. Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control unit or other SCB units (see SCB Instruction Manual).

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Operating the unit

Operating the unit for the first time

The AUTOCON® II 400 is delivered as follows:

Language: English

Neutral electrode: universal

Maximum activation time: 99 sec.

Predefined procedures

For every unit 100 user procedures can be saved. The memory locations 0 to 99 are reserved for this purpose. Every unit is equipped with several predefined procedures which cannot be modified or deleted. They occupy memory locations ab 100.

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Default values of HF modes:

Operating mode designation	Default effect value	Default output value	Default Setting
Cut Off	-	-	
TOP-Cut	1	50	
POWER-Cut	1	50	
Bipolar-Cut	1	40	
C-Cut	1	140	
Lap-C-Cut	1	50	
Saline-C-Cut	6	-	
Saline-Time-C-Cut	6	-	time 0.2 sec.
Gastro-Cut	3	120	ECUT-SPEED = 8 ECUT-IMPULSE = 2
Papillo-Cut	2	100	ECUT-SPEED = 8 ECUT-IMPULSE = 2
Coag Off	-	-	
Standard Coag	1	40	
Forced Coag	1	40	
Spray Coag	1	40	
Bipolar Soft Coag	1	25	
Bipolar Soft with Auto-Stop	1	25	
Bi-Vascular-Safe	1	30	
Saline-Coag	1	-	
Saline-Time-Coag	6	-	time 0.3 Sec.

Название рабочего режима	Стандарт. значение эффекта	Стандартное значение ограничения мощности	Заводская настройка
Cut Off	-	-	
TOP-Cut	1	50	
POWER-Cut	1	50	
Bipolar-Cut	1	40	
C-Cut	1	140	
Lap-C-Cut	1	50	
Saline-C-Cut	6	-	
Saline-Time-C-Cut	6	-	0,2 с.
Gastro-Cut	3	120	ECUT-SPEED = 8 ECUT-IMPULSE = 2
Papillo-Cut	2	100	ECUT-SPEED = 8 ECUT-IMPULSE = 2
Coag Off	-	-	
Standard Coag	1	40	
Forced Coag	1	40	
Spray Coag	1	40	
Bipolar Soft Coag	1	25	
Bipolar Soft con Auto -Stop	1	25	
Bi-Vascular -Safe	1	30	
Saline-Coag	1	-	
Saline-Time-Coag	6	-	0,3 с.

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Switch on the power switch ①.

When the unit is switched on, the logo screen is first displayed for approximately 3 seconds, then a white screen appears for approx. 2 seconds.

During power-up the unit performs a self-test. Here it checks important system components such as the HF generator, output sockets, touch-screen, internal memory components, etc. and indicates any malfunction via an error message. The two signal paths for alarms, display and the audio system must be checked by the user whenever the unit is powered up.

Note: During power-up of the unit, the KARL STORZ logo must firstly be displayed followed by the white screen. In the event of any error the unit must not be used.

Warning: When powering up the unit, pay attention to the acoustic availability signal.

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
The power-up message is followed by the “Procedure List” screen. This lists the predefined user procedures and the procedures defined by the user himself. In routine operation the overview of procedures is called up here. If button ① is activated, the user will see the screen “Device Setup: Select Option” (when the unit is operated for the first time, English is preset).

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
Selection of language and user settings

The button “User Settings” ②4 can be used to access the selection list of the user settings.

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The arrow keys ② next to the field “Language” can be used to select the language. With this menu the user can also select the brightness of the touch-screen, the volume of the alarm system and the activation signals, the keyclick volume and the date format. If the user activates button ③, the properties selected will be permanently saved.

If the user activates button ④, he will return to the “User Settings” screen without saving the values. If the user continues to activate the same button, he will return to the “Procedure List”.

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
Selecting and modifying a predefined procedure

One of the predefined procedures can be selected from the menu “Procedure List” by activating the relevant assignment buttons.

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Note: With the procedure button ⑱ the user will return to the procedure list.

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An overview of sockets will appear with the relevant parameters.

To modify the operating parameter, activate the button of the relevant operating mode, “Cut” (yellow) or “Coag” (blue).

Note: The buttons ⑤, ⑥ and ⑦ allow the user to switch between the individual sockets.

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
The parameters can now be adjusted with the Plus/Minus buttons. The newly selected values are active but not permanently saved. They can be permanently saved by following the instructions specified under the item "Saving a personal self-defined procedure" (p. 54).

Note: With the button ④ the user will return to the "overview of sockets".

With the ⑯ mode button the user can access mode selection.

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
Here the user can make his selection as required for the relevant socket. The different sockets offer their own selection for the relevant modes "Cutting" and "Coagulation". Activating the button will switch the unit to Active mode. This is indicated by the green indicator in the left-hand corner of the button. The unit automatically switches to the overview of sockets.

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
In this overview, the buttons ⑤ and ⑦ can be used to select the output parameters of the other sockets. The socket activated in each case is indicated by a green indicator in the top left corner of the button. At the same time, the pictogram showing zugewiesene Aktivierungsart ⑨/⑩ is also lit up. The neutral electrode button ⑧ can be used to access the Neutral Electrode menu.

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Here the contact quality of any neutral electrode applied is given in ohms. An acceptable value is shown in green, an unacceptable value in red.

The button ④ is used to navigate back to the overview of sockets, and from there back to the overview of procedures.

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
Saving a personal self-defined procedure

If parameters have been amended in a “predefined procedure” as instructed under the steps “Selecting and modifying a predefined procedure” (p. 52), they can be permanently saved for subsequent applications. The button ⑪ should be pressed for this purpose.

How to overwrite and delete a procedure is described from p. 65 onwards.

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The button “Save Procedure As” ⑮ can be used to access the menu to save the procedure.

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The Editor can be accessed by directly selecting the memory location or navigating with the arrow keys and then selecting the memory location. The number of the memory location is shown on the right in the field for the button.

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
Note: If upper-case letters have been activated, the green indicator on the shift key ⑰ lights up.

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Saving a level “B” with a self-defined procedure

With the optional triple-pedal footswitch **20013831** it is possible to switch backwards and forwards between two procedures of an application as required. Here a self-defined procedure as described under “Selecting and modifying a predefined procedure” (p. 52) should be created. **Diese Prozedur ist durch Auswahl neuer Parameter und Fußschalterfunktionen entsprechend den Bedürfnissen zu verändern.** When saving as described under “Saving a personal self-defined procedure”, the button “Create Level B of Procedure” should be pressed in the menu to save a procedure. This allows the user to access the Editor for level “B” of a procedure.

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The two procedures will now appear one under the other in the procedure list.

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
In the overview of sockets of a procedure the relevant socket is assigned a footswitch activation. By adjusting or modifying a procedure in line with current requirements, the relevant footswitches, YELLOW for cutting and BLUE for coagulation, which are shown in the menu by colored pictograms next to the socket buttons, can be reassigned to another socket or blocked. To do so press the relevant button ⑤, ⑥ or ⑦ next to the socket. The button which can be activated has a light gray background and is assigned to the relevant socket. It is possible to connect two footswitches to the unit: a single-pedal footswitch and either a dual or triple-pedal footswitch. The single-pedal switch and thus also its pictogram only have one blue pedal. This means that it is only possible to activate the coagulation function.

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The footswitch activation selected, here blue for coagulation, is indicated by the green indicator ⑨. After selection the screen view automatically switches back to the overview of sockets, where footswitch activation can be assigned to additional sockets or blocked. With barred footswitch activation, HF energy can be released in monopolar modes via the buttons on the instrument (hand switch). For bipolar applications automatic activation can be selected. This is time-dependent and can be specified in the menu "Device Setup".

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Autostart function

The Autostart function permits in the coagulation modus of the bipolar socket, automatic activation of the coagulation current as soon as both branches of the coagulation electrode touch the tissue and if there is a defined impedance area. The Autostart function 1 allows coagulation to be activated deferred from between 0.1 and 1.9 seconds in steps of 0.1 second. The Autostart function 2 provides for a time delay of between 2 and 10.0 seconds in steps of 0.5 seconds.

Autostart 1 is factory-set to a value of 1.0 seconds, and Autostart 2 to a value of 2.0 seconds. The user is free to reselect the time-to-autostart delays at his own discretion.

Note: The Autostart function can be used for the operation modes

- Bipolar Soft Coag
- Bipolar Soft Coag with Auto-Stop.

To select the autostart function or end delay time, a bipolar HF socket must be selected on the touch screen.

Following selection this is confirmed by the green display (5).

The field for selection of the activation type is now shown next to the pictographic HF socket.

The activation times 1 and 2 are now available for selection in the selection list next to footpedal activation.

The permitted time selected under "Application Settings" (see p. 58) is then shown.

Warning: The HF generator is now automatically activated when the tips of the forceps come into contact with tissue!

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- Bipolar Soft Coag
- Bipolar Soft Coag con Auto-Stop.

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Time setting for bipolar Autostart activation

To set the time for bipolar Autostart activation, the button ① "Setup Device" should be activated from the "Procedure List" screen.

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The button "Application Settings" ⑱ can be used to switch to the overview in which application settings are made.

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
The required times can be set with the arrow keys next to the fields Autostart 1 and Autostart 2. After selection, button ③ is used to save the parameters, and the system automatically returns to the menu "Device Setup: Select Option". With button ④ the user can move back without saving changes.

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Selecting the neutral electrode

The type of neutral electrode is selected for all monopolar procedures of the AUTOCON® II 400. It can be adjusted as follows. From the start screen "Procedure List" activate the button ① "Setup Device".

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The button “Application Settings” ⑱ can be used to switch to the overview in which application settings are made.

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The button “NE Type” can be used to switch to the selection list for the neutral electrode. The button “NE Type” indicates the type of neutral electrode selected.

The required mode can be selected with the arrow keys next to the field NE Type. The settings single-faced, double-faced, dynamic and ‘any’ can be selected for the connection of single-faced and double-faced neutral electrodes. The measuring range is from 0 - 199 ohms (lower and upper ends of bar). Greater resistances are shown in the numerical display by “---”. Activation of HF energy (p. 62) on this level will take the user to level 2 (procedure information) as not all parameters are displayed. HF energy will only be output with renewed activation on this level.

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Setting the maximum activation time

The maximum activation time can be adjusted according to the user's requirements. From the start screen "Procedure List" activate the button ① 'Setup Device'. This setting is protected with a user password.

Warning: Any change to the maximum activation time must be notified to all users.

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The button ⑱ "Application Settings" can be used to switch to the overview in which application settings are made.

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The user can use the arrow keys ⑲ next to the field "Max. Activation Time" to set the required time. When the arrow keys ⑲ are activated, the menu for entering the user password will firstly appear.

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Note: On delivery the user password is "12345". It is not identical with the service password and can be changed as described in the section "Changing the user password" (p. 65).

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Installation and operating instructions

Entry of the password is secret, i.e. only asterisks appear in the display in place of the characters entered. After entering the valid password the user can confirm using button ③ and return to the “Application Settings” menu where the buttons for selection of the maximum activation time have been enabled. After specifying the activation time this value can be permanently saved with button ③. To make any change to the activation time the valid password must be reentered.

Connecting instruments

Depending on the type of unit, the monopolar and/or bipolar instruments can be connected to the sockets on the front using the appropriate cable.

Connecting the footswitch

At the rear of the unit the user can connect the following as required:

- a single-pedal footswitch,
- a dual-pedal footswitch,
- a triple-pedal footswitch or
- a combination of a dual-pedal footswitch and a single-pedal footswitch and
- a combination of a triple-pedal footswitch and a single-pedal footswitch.

The ‘Assigning **Aktivierungsart** to a socket’ is described on p. 56.

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Error messages

The AUTOCON® II 400 is equipped with an error control system. User errors and unit errors are flashed on the screen (13). The error codes are listed in the table on page 81.

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Further setting options:

Setting date and time

The AUTOCON® II 400 is supplied with the date and time set to 'Central European Time' CET. From the start screen 'Procedure List' activate the button ① Setup Device.

Note: The ticking clock shows that the screen is being updated. If the clock has stopped, the screen is frozen, and this error must be rectified before using the HF unit again.

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With the button ⑳ “Date and Time” the user can switch to the menu used to set the date and time.

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The settings can be made using the arrow keys under the fields. After selection, button ③ is used to save the values, and the system automatically returns to the menu “Device Setup: Select Option”.

Note: The day of the week is calculated using the calendar. Selection is not necessary here.

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
The format for display of the date can also be selected. With the button “User Settings” ⑳ the user can access the user settings menu.

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If the user activates button ③, the properties selected will be permanently stored. If the user activates button ④, he will return to the screen “User Settings”. If the user activates this button again, he will return to the “Procedure List”.

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
Changing the user password

If button ① is pressed from the menu "Procedure List" the user can access the screen "Device Setup: Select Option".

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The button "Change User Password" ②5 can be used to access the menu "Enter / Edit Password".

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Entry of the password (the password has been factory-set to "12345") is secret for security reasons - ***** appears in the display in place of the characters entered. The entry of a valid password takes the user to another entry menu, which only differs from the previous menu by having two entry lines. The text of the header changes to 'Input New Password'.

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Note: The user password is not identical with the service password.

After the new password has been entered in line 1, pressing the button ⑫ will move the cursor to the second input line.

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The text in the header changes to “Confirm New Password”. After confirming the password, pressing button ③ will take the user back to “Device Setup: Select Option” if the entry is valid, and the new password will be accepted. From here you can return to the menu “Procedure List” with button ④. Activation of the button ③ with an invalid password will result in the message “Invalid password”. The content of the two text lines is then deleted and the cursor is again positioned in the left-hand input line. Pressing the Back button will take the user back to the Set-Up selection list, with the password remaining unchanged.

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
Overwriting a procedure

Here there is the option of overwriting a procedure if the parameters of a procedure have not proved to be effective and are to be permanently amended. The procedure to be overwritten is selected from the procedure list by pressing the appropriate button. This accesses the overview of procedures.

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Parameters are selected as described under “Selecting and modifying a predefined procedure” (p. 52). From the overview of procedures, button ⑪ can be used to access the menu “Edit Procedure”.

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The button “Overwrite Procedure” is used to determine the procedure displayed in the header to be overwritten. The procedure is overwritten using button ③, and the user is then taken to the procedure list.

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
Deleting a procedure

The procedure to be deleted is selected from the procedure list by pressing the appropriate button. This accesses the overview of procedures.

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From the overview of procedures, button ⑪ can be used to access the menu "Speicheroptionen".

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The button "Delete Procedure" ⑳ is used to determine the procedure displayed in the header to be deleted. The procedure is deleted using button ③, and the user is then taken to the procedure list.

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Display of system events

The system events are listed in the logbook. From the "Procedure List" menu, button ① can be used to access the menu "Device Setup: Select Option".

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
The button “Show Log” 28 is used to access the logbook in which the system events are stored.

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
The system events can be viewed using the arrow keys on the right of the Fig. With button 4 the user can return to the menu “Device Setup Select: Option”.

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From here you can return to the menu “Procedure List” using button 4.

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Display of device information

The device information states the serial number and software versions. From the menu “Procedure List”, button 1 can be activated to access the menu “Device Setup Select Option”.

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
The button “Device Configuration” 29 is used to access the device information where the configuration data is stored.

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With button ④ the user can return to the menu "Device Setup: Select Option".

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From here the user can return to the menu "Procedure List" using button ④.

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Service Programs

The button "Service Programs" ⑩ can be used to access the password-protected service area. For Technical Service the service menus allow certain system information to be displayed: to reset the date of the next safety check and allocate a new service password. The service password is not identical with the user password. Further details can be found in the Service Manual for this device.

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In this menu it is possible to select the activation times (0.1 s - 1 s, in 100 ms steps) for the modes Saline-Time-Cut and Saline-Time-Coag using the UP/DOWN buttons. Activation is stopped according to the time selected regardless of whether the footswitch or fingerswitch is activated again.

In addition, the number of cutting pulses (1 - 5) and the cutting speed from slow to fast (1 - 20) can be selected separately for Gastro- and Papillo-Cut.

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Fuse replacement

Switch off the unit and remove the power cord from the electrical outlet.

Only use fuses suitable for the voltage range stated on the manufacturer's identification plate.

Caution: Use only fuses of the correct rating. Replace defective fuses (see table) and reinsert line fuse holder/voltage selection module in the unit.

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Release cover from the line fuse holder/voltage selection module with a screwdriver or other suitable tool. Remove the line fuse holder/voltage selection module.

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Caution: Only use fuses of the correct rating.

Replace defective fuses (see table) and reinsert line fuse holder/voltage selection module in the unit.

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Reconnect the power cord. Test the unit for proper operation.

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Cleaning and sterilization

AUTOCON® II 400

Warning: Always pull out power plug before cleaning!

Caution: Avoid allowing fluids to enter the unit.

The exterior of the unit and footswitch can be cleaned with a conventional surface disinfectant.

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Manual/machine cleaning and disinfection

Note: Use a cleaning/disinfecting agent approved by KARL STORZ for cleaning/disinfection (see Appendix, page 109). You will also find the current list of approved products on our Web site (www.karlstorz.com). Further information can be found in the general instructions for the preparation of instruments.

Caution: Follow the solution manufacturer's instructions, paying close attention to proper dilution and recommended minimum exposure times. Prolonged exposure may result in corrosion or material damage.

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Sterilization of instruments

Note: Before sterilization, the components must be cleaned:

- mechanically,
- using ultrasound or a
- washing machine (up to 93°C)

Then pack the components into appropriate containers.

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Warning: Sterilization is only possible on clean surfaces.

Warning: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.


Warning: Any deviations from the recommended parameters for sterilization should be validated by the user

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Steam sterilization

Sterilization trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

Load the sterilizer according to its instruction manual.

On completion of the steam sterilization cycle, all sterilized parts should be allowed to cool down slowly.

Warning: Risk of burns! Parts are hot after steam sterilization. Allow them to cool down!

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Validated steam sterilization methods


KARL STORZ has validated the following steam sterilization methods:

- Prevacuum
- Fractionated prevacuum
- Gravity displacement

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Vacuum/prevacuum method

The prevacuum method consists of four phases.

- **Conditioning phase**
In the conditioning phase a vacuum is created in the sterilizing chamber; the instruments are then heated by injecting steam.
- **Sterilization phase**
Sterilization takes place at a temperature of 134°C and an operating pressure (p_b) of 2 bar over a period of 5-8 minutes.
- **Evacuation phase**
The steam is released.
- **Drying phase**
Drying takes place under renewed application of a vacuum over a period of approx. 5-20 minutes.

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
Fractionated prevacuum method

The fractionated prevacuum procedure also consists of the four previously mentioned phases. In the conditioning phase, a vacuum is generated up to four times. Steam is then blown in and sterilization takes place at 121°C (250 to 252°F) or 134°C (270 to 272°F) in 20 or 5 minutes respectively (minimum times). Before the final drying phase, the steam is removed from the chamber during the evacuation phase. Drying takes place under vacuum.

Validated parameters for the fractionated prevacuum method:


Temperature:	270 to 272 °F; (134 °C ± 3 °C)
Pressure (p_b):	2 bar
Exposure time:	5 minutes

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Gravity displacement method

Warning: This method is only suitable for smooth and accessible surfaces, not for instruments with channels, gaps or stop-cocks.

The gravity displacement method is also made up of 4 phases:

- **Conditioning phase**
Steam is pressed into the chamber and heats the instruments.
- **Sterilization phase**
Sterilization takes place at a temperature of 121°C and an operating pressure (pb) of 1 bar over a period of 30 - 45 minutes.
Sterilization may also take place at 134°C and 2 bar; the time required is 5 - 10 minutes.
- **Evacuation phase**
Steam is released through the outlet.
- **Drying phase**
Drying takes place under ambient air pressure by heating through the heating jacket of the autoclave.

Validated parameters for the gravity displacement method:

Temperature:	270° to 272°F, 121°C
Pressure (pb):	1 bar
Required time:	30 - 45 minutes

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Maintenance

Cleaning and sterilization

* Gas mixture composed of 10% ethylene oxide and
90% chlorotetrafluoroethane (HCFC-124)
(percent by weight)
Airing can be accomplished at a temperature of
50–55°C in 12 hours. This increases to 24 hours at
room temperature (20°C).

Maximum permissible levels (in accordance with 21
CFR 812.100; US Federal Register, Vol. 43, No. 122,
June 23, 1978) for residues following EtO steriliza-
tion are as follows:

Ethylene oxide: 25 ppm* / 250 ppm**

Ethylene chlorohydrine: 25 ppm* / 250 ppm**

Ethylene glycol: 250 ppm* / 5000 ppm**

** for instruments in contact with skin and
mucous membrane,

* for instruments in contact with blood.

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Gruendemann, B.J. and Meeker, M.H. Alexander's
Care of the Patient in Surgery, 7th edition. The C.V.
Mosby Company, St. Louis, Mo. 1983.
Association for the Advancement of Medical
Instrumentation. Designing, Testing and Labeling
Reusable Medical Devices for Reprocessing in
Health Care Facilities: A Guide for Device
Manufacturers. AAMI TIR No. 12-1994.
Descoteaux, J-G, Poulin, E.C., Julein, M. and
Guidoin, R. Residual Organic Debris on Processed
Surgical Instruments. AORN Journal,
Vol. 62, No. 1, pp. 23-29, 1995.

Gas sterilization with formaldehyde (FO)

Caution: During formaldehyde steriliza-
tion the temperature must not exceed
65°C.

When using formaldehyde the instruments can be
used immediately after the sterilization cycle has
been completed. An additional airing time is not
necessary.

Note: Observe national laws and regulations for the
gas sterilization methods.

References for Cleaning and Sterilization

Gruendemann, B.J. and Meeker, M.H. Alexander's
Care of the Patient in Surgery, 7th edition. The C.V.
Mosby Company, St. Louis, Mo. 1983.
Association for the Advancement of Medical
Instrumentation. Designing, Testing and Labeling
Reusable Medical Devices for Reprocessing in
Health Care Facilities: A Guide for Device
Manufacturers. AAMI TIR No. 12-1994.
Descoteaux, J-G, Poulin, E.C., Julein, M. And
Guidoin, R. Residual Organic Debris on Processed
Surgical Instruments. AORN Journal Vol. 62, No. 1,
pp. 23 - 29, 1995.

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
Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

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Safety checks and inspections

In order to prevent accidents due to aging, wear or failure of medical equipment, regular safety checks are recommended.

The high frequency surgical unit must be inspected at least once a year. Please observe the accident prevention regulations and testing intervals applicable in the various countries.

12 months after the last safety check and inspection, a status message is displayed when the unit is switched on. This appears briefly and serves as a reminder.

Warning: Servicing may only be performed by specialists authorized by KARL STORZ.

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If any defects or shortcomings which could endanger patients, operating personnel or others are ascertained during these safety checks, the unit must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

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In accordance with the Medical Device Directive (MDD) this instruction manual contains information on the type, extent and frequency of safety checks.

When carrying out safety checks, the following items must be tested:

- Visual inspection of the unit and accessories for proper condition
- Checks of protective grounding and insulation, and for leakage currents
- Testing for proper operation of all visual and acoustic signals
- Testing of the monitoring devices
- Testing of the function modes
- Testing of the automatic start mode
- Checking the activation of the HF power using the dual-pedal footswitch
- Testing the HF power in the various operating modes
- Measurement of power consumption when operating at no-load and under load

Note: For detailed information on the range and implementation of the safety checks please refer to the Service Manual.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

For fiberscopes and equipment, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen/Germany

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen/Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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AUTOCON® II 400 Error Messages

Error code	Symptom	Troubleshooting	Priority
A-0B:	NE contact: Resistance too low	Check neutral electrode and reposition if necessary	Medium
A-0C:	NE contact: Resistance too high	Check neutral electrode and reposition if necessary	Medium
A-10:	End activation (AutoStop reached)	End activation of HF energy	Low
A-11:	No activation signal		High
A-12:	Activation during power up not permitted	Check footswitch, handle and cable	Medium
A-17:	Simultaneous activation not permitted	Activate either via footswitch or handle	Medium
A-19:	Line voltage fault present	Check power cord, power outlet and power supply	Low
A-1C:	Max. activation time exceeded	Cancel activation	Low
A-1E:	Touch screen pressed during power up	Touch screen must not be pressed during power-up	Low
A-1F:	NE symmetry critical (cancel activation)	Check neutral electrode and reposition if necessary	High
A-22:	AutoStop reached	Break contact with tissue	Low
A-85:	Dual-pedal footswitch disconnected	Reconnect dual-pedal footswitch and check cable	Low
A-89:	Single-pedal footswitch disconnected	Reconnect single-pedal footswitch and check cable	Low
A-94:	Triple-pedal footswitch disconnected	Reconnect triple-pedal footswitch and check cable	Low
A-95:	Instrument an MF-Buchse erkannt	New MF instrument at socket	Low
A-99:	Activation type not available	Activation type not available with instrument connected	Low
A-9F:	Triple-pedal instrument disconnected, socket n	MF instrument disconnected	Low
A-A3:	Footswitch not assigned to a socket	Assign footswitch to a socket (p. 56)	Low
A-A4:	2 footswitches of same type connected	Either connect one footswitch or two of a different type	Low
A-A9:	First select procedure before activation	Select the required procedure (p. 52)	Low
A-AA	No mode assigned	Assign a mode to the procedure (p. 30)	Low
A-AB	Instrument detection defective	Check instrument and cable	Medium
A-AC	Autostart on tissue contact not possible	Release tissue contact, Instrumentbranchen ggf. trennen, Autostart zuordnen	Medium
A-B0	NE symmetry not optimal	Check neutral electrode and reposition if necessary	Low
A-B7	Autostart only permitted up to 50W output	The max. output with Autostart should only be 50W	Low
A-BB	SC due	Arrange date for safety check	Medium
A-C0	Activation type for MF-instrument missing	Specify activation type	Low
A-FE	Line voltage fault present	Check power cord, power outlet and power supply	Low
C-75	Error with instrument detection	Check instrument and instrument cable and reconnect	Low
C-84	Short-circuit between AE and NE	No impedance modification of tissue during coagulation	Low
I-84	Defective instrument detection	Check instrument and instrument cable and reconnect	Low
I-8B	Impermissible instrument detected	Check connected instrument and replace if necessary	Low
I-8D	Impermissible instrument detected	Check connected instrument and replace if necessary	Low
I-8E	Impermissible instrument detected	Check connected instrument and replace if necessary	Low
all others		Inform Service	Low

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POWER Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	740 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	10 W to 300 W in 1 W steps
Max. power output at rated load resistance	300 W ± 20%
TOP Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	950 Vp (with electric arc)
Number of effects	8
Constancy of effects	automatic control of arc intensity
HF power limitation	10 W to 300 W in 1 W steps
Max. power output at rated load resistance	300 W ± 20%
LAP-C-Cut	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	
Effect 1-4:	3 (at RL = 500 ohms)
Effect 5+6:	3.2 (at RL = 500 ohms)
Effect 7+8:	3.8 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	1450 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	10 W to 200 W in 1 W steps
Max. power output at rated load resistance	200 W ± 20%

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PAPILLO Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Max. HF peak voltage	800 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
Max. power output	200 W ± 20%
GASTRO Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Max. HF peak voltage	880 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
Max. power output	200 W ± 20%
BIPOLAR Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	740 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
Max. power output	100 W ± 20%
High frequency power limitation	1 W to 100 W in 1 W steps
C-Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	
Effect 1-4:	3.0 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	1450 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
Max. power output	200 W ± 20%
High frequency power limitation	10 W to 200 W in 1 W steps

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SALINE-C-Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	640 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
Max. power output	370 W +8%/-20%
SALINE-TIME-C-Cut	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	640 Vp
Number of effects	8 (time 0.1 - 1sec)
Constancy of effects	automatic control of HF peak voltage
Max. power output	370 W +8%/-20%
STANDARD Coag	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	50 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	5 W to 200 W in 1 W steps
Max. power output at rated load resistance	200 W ± 20%

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FORCED Coag	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	6.0 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	1800 Vp
Number of effects	4
Constancy of effects	automatic control of HF peak voltage
HF power limitation	5 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W ± 20%
SPRAY Coag	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	7.4 (at RL = 500 ohms)
Rated load resistance	500 ohms
Max. HF peak voltage	4300 Vp
Number of effects	2
Constancy of effects	limitation of HF peak voltage
HF power limitation	5 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W ± 20%
BIPOLAR SOFT Coag / BIPOLAR SOFT with AUTOSTOP	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency 350 kHz	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	75 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	1 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W ± 20%

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BI-VASCULAR-Safe	
HF voltage form	unmodulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	200 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
HF power limitation	1 W to 120 W in 1 W steps
Max. power output at rated load resistance	120 W ± 20%
SALINE Coag	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	75 ohms
Max. HF peak voltage	190 Vp
Number of effects	8
Constancy of effects	automatic control of HF peak voltage
Max. power output at rated load resistance	200 W ± 20%

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SALINE-TIME-Coag	
HF voltage form	pulse-modulated sinusoidal AC voltage
Rated frequency	350 kHz (at RL = 500 ohms) ± 10%
Crest factor	1.4 (at RL = 500 ohms)
Rated load resistance	75 ohms
Max. HF peak voltage	190 Vp
Number of effects	8 (time 0.1 - 1sec)
Constancy of effects	automatic control of HF peak voltage
Max. power output at rated load resistance	200 W ± 20%

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






























Technical description**Technical data****Standard compliance** (for 205352 20-11x)**According to IEC 60601-1; 60601-2-2:**

- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type CF 

Directive compliance (for 205352 20-11x)

This medical device belongs to Class IIb
This medical device bears the CE mark in accor-
dance with the Medical Device Directive (MDD)
93/42/EEC.
A code number after the CE mark indicates the
responsible notified body.

Technical documentation

On request, the manufacturer will provide those
circuit diagrams, itemized parts listings, descrip-
tions, sets of adjustment instructions, and other
items of available documentation to suitably
qualified user personnel duly authorized by the
manufacturer for their use in repairing those
components of the unit that have been designat-
ed by their respective manufacturers as
repairable.

Supply of such technical documentation relating
to the unit shall not be construed as constituting
manufacturer's authorization of user's personnel,
regardless of their level of technical training, to
open or repair the unit.

Explicitly exempted herefrom are those mainte-
nance and repair operations described in this
manual.

Note: We reserve the right to make engineering
modifications in the interest of promoting techno-
logical progress and generating performance
improvements without obligation on our part to
submit prior notice thereof.

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
**Spare parts,
recommended accessories**

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Chemicals for the preparation of telescopes and endoscopic instrument sets
Manual cleaning and disinfection/ultrasound

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, motors and their handles, IMPERATOR drill handles, prismatic light deflectors. The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to proper dilution and recommended immersion times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept forte	2, 3
	Sekusept Plus	1, 2, 3
	Sekusept Pulver Classic	1, 2, 3
	Sekusept Pulver Classic + Aktivator	1, 2, 3
Esteer GmbH	Ultrasdesmit Konz.	2, 3
Fresenius AG	Afid	1, 2
	Afid plus (neu)	1, 2
	Sporcid FF	1, 2
Indeba Industria E Comercio LTDA	Letahdeido	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehyd 2000	2
Lysoform Dr. Rosemann GmbH	Aldasan 2000	1, 2
	Almyrol	1, 2
	Desoform	1, 2
	Lysoformin 3000	1, 2
Merz+Co.	Edisonite Super	3
	Mucadont-IS	2, 3
	Mucadont-Zymaktiv	1, 3
	Mucoicid-T	1, 2, 3
Orochemie	A 10	2
	ID 210	2
Promagent AB	Wavacide	2
Schülke & Mayr GmbH	Gigasept	2, 3
	Gigasept FF	2, 3
	Lysetol AF	1, 2, 3
	Lysetol V	2
Schuhmacher, Dr. GmbH	Descoton forte	2
	Perfektan Endo	2
Steris	EnzyCare 2	1, 3
	neodisher LM 2	1
Weigert, Dr. GmbH & Co.	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher medizym	1, 3
Whiteley Industries PTY.Ltd	Aidal Plus	2

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination shaft, **not** the housing/handle part.

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths

Note: You will also find the current list of approved products on our Web site (www.karlstorz.de or www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid

Note: You will also find the current list of approved products on our Web site (www.karlstorz.de or www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine cleaning and disinfection

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, motors and their handles, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schuhmacher, Dr. GmbH	Thermoton Endo	2
Steris	Instru-Klenz	1, 3
	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher mediclean	1
	neodisher mediclean forte	1
	neodisher mediklar	5
	neodisher medizym	1
	neodisher Septo DA	2
	neodisher Septo DN	2
	Weigert, Dr. GmbH & Co.	

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this Appendix, as well as the instruction manual, during installation and commissioning.

The AUTOCON® II 400 model 20 5352 20-11x corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The AUTOCON® II 400 model 20 5352 20-11x is a Group 2 unit (as per CISPR 11).

Group 2 contains all the 'equipment and systems which generate **electromagnetic energy** to perform their intended function. **Adjacent electronic devices may be affected**'.

NOTE: The tables and guidelines that are included in this Appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The AUTOCON® II 400 model 20 5352 20-11x should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the AUTOCON® II 400 model 20 5352 20-11x. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200				
Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
Handle BOWA 214-03	No		No	Handle
205300 76	No		No	Neutral electrode
Handle MTP	No		No	Handle
Erbe 20489-100	Yes	4.5 m	No	Double pedal footswitch
Erbe 20188-100	Yes	4.5 m	No	One pedal footswitch
200901 70	No	5.0 m	No	ME Cord
205300 70	No	3.0 m	No	NE Cord
200901 70	Yes	1.0 m	No	SCB Cord
PE	No	3.0 m	No	Potential equalization
Power Cord	No	3.0 m	No	Power Supply

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Table 201		
Guidance and manufacturer's declaration – electromagnetic emissions		
<p>AUTOCON® II 400 model 20 5352 20-11x is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® II 400 model 20 5352 20-11x should ensure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	<p>The AUTOCON® II 400 model 20 5352 20-11x is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p>
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000 -3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202

Guidance and manufacturer's declaration – electromagnetic immunity

AUTOCON® II 400 model 20 5352 20-11x is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® II 400 model 20 5352 20-11x should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 204

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

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Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V_{rms} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V_{rms}</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the AUTOCON® II 400 model 20 5352 20-11x, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = [3.5/3]\sqrt{P}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
AUTOCON® II 400 model 20 5352 20-11x**

The AUTOCON® II 400 model 20 5352 20-11x is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

To be completed by the supplier/importer:

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To be filled out by instrument owner:
Return address/company stamp: _____

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KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Rm 1119-20, Dongshan Plaza
69 Xianlie Road (Middle),
Dongshan District,
Guangzhou, PRC., 510095 China
Phone: +86 10 8732-1281
+86 10 8732-1282
Fax: +86 10 8732-1286

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966
Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

*Sales for Industrial Endoscopy

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com