Better Vision Better Results

WIRELESS

With a dedicated wireless link and a 3.5 hour Lithium battery, DE1250 provides unparalleled flexibility in performing endoscopic observations and capturing high precision videos and images.

100x MAGNIFICATION

Provides up to 35x native optical magnification and up to 100x digital magnification

LIGHTWEIGHT

Based on a lightweight but rugged construction, the DE1250 is easy to use without hindering mobility during endoscopy.

MULTI-LAYER GLASS LENSES

Equipped with multi-layer glass lenses rather than plastic alternatives, thus delivering superb focal depth and crystal clear image quality.



- Completely wireless operation
- Compatible with rigid endoscopes
- Captures snapshots & videos
- Ideal for ENT, teaching and veterinary applications
- Magnifies objects up to 35x optically and up to 100x digitally

Firefly DE1250 is the industry's first wireless compact digital Endoscope Camera with image and video capture capabilities. It provides unprecedented viewing convenience with popular rigid and flexible endoscopes. Besides endoscopic examinations, this powerful tool can also be used for teaching applications, for EMR (Electronic Medical Records) as well as for patient education. Unlike older endoscope cameras, the DE1250 delivers these breakthrough capabilities - and much more - at a highly economical cost.



Copyright © 2014 by Firefly Global, All rights reserved. Rev 1.5 www.fireflyglobal.com

Better Vision Better Results

ADVANCED IMAGING **TECHNOLOGY**

Managing incoming light is a tremendous challenge in digital sensor design. Reflections and lighting leaks can easily pollute the sensor, resulting in overexposure and washed-out images. Firefly effectively reduces overexposure by running highly efficient image processing algorithms and by employing an advanced lens design.

Sophisticated image processing algorithms are ineffective without strong depth perception and a clear image. Firefly DE1250 is equipped with multi-layer glass lenses rather than basic alternatives, thus delivering superb focal depth and crystal clear images.

The DE1250 also leverages its advanced CPU to the fullest extent: Videos and images are transferred to the PC in realtime, up to 30 frames-persecond for a live viewing experience.

DE1250

Technical Details

Magnification Native Optical: 15x - 35x Digital: 15x - 100x Sensor Resolution 720 x 480 Lens Assembly Dual Lenses, 3-Layer Glass F = 22mmVideo Format: YUY2 Frame rate: 30 FPS Video/Image Properties Hue, Saturation Color: Exposure: Brightness, Contrast Image: Sharpness, Gamma Image / Video Files JPG / BMP / AVI Wireless Link 4 Channels Built-in 850mAH Battery **Dimensions** 11.5cm Length 3.6cm Rear Diameter 4.6cm Coupler Diameter Interface USB 2.0 AC Battery Charger Accessories **Operating Systems** Win 8, Win 7, Vista and XP Mac OS-X 10.4 or higher Software Scalable Window, Zoom, Freeze, Resolution, Rotate, Flip

INCLUDES WIRELESS RECEIVER

Region of Interest (ROI)

1 Year Limited Hardware Warranty



Specifications are subject to change in any matter and at any time without notice. The Firefly word and the Firefly logo are trademarks or trade names of Firefly Global and its affiliates in the United States and/or other countries. All rights in such names, marks or logos are reserved by Firefly Global.

Copyright © 2014 by Firefly Global, All rights reserved. Rev 1.5

www.fireflyglobal.com

Warranty

Firefly



SOPHISTICATED SOFTWARE

DE1250 is controlled directly with the FireflyPro professional image processing software which is bundled with the product. This state-of-the art software enables users to capture, store, recall, view, manipulate and measure images and videos in real time. Its intuitive interface empowers users to get right to work anytime and anywhere.



OPTIONAL TV KIT

Using the optional ES160 TV Kit, a standard LCD/LED TV can also be used to display images and videos being streamed from the DE1250.